

AUDIBEL  iQ  
CROS SYSTEM

---

Hearing solutions for  
single-sided hearing loss



  
AUDIBEL™

# Say hello to a solution for single-sided hearing loss

Audibel uses cutting-edge technology to deliver true listening enjoyment.

With our new line of **Audibel A4™ iQ** hearing aids, Audibel™ has a solution for those with single-sided hearing loss.



AUDIBEL **A4** iQ

*Audibel offers the first rechargeable hearing aid for single-sided hearing loss.*

Now, if you have single-sided hearing loss you'll be able to hear the world around you in full, pristine sound.

## *The A4 iQ CROS System:*

1. Allows you to hear from both sides — with devices so small and discreet they're nearly invisible
2. Improves listening in noise
3. Reduces background noise



# Say hello to customised solutions



*Also available with a*



## CROS solution

If you're unable to hear in one ear and have normal hearing in the other ear, you may benefit from a CROS solution. Finally, you'll be able to hear sounds in your weaker-hearing ear by way of your healthy ear.

### How it works

To help you better engage in all types of listening environments, our CROS solution has a microphone that picks up sounds and voices from your weaker-hearing ear, then transmits them to a hearing aid receiver that's fitted on your good ear. No matter your position at the dinner table, or if sounds originate from your weaker-hearing side, you won't miss important sounds around you.




## BiCROS solution

A BiCROS solution is also available and can benefit those with little to no hearing in one of their ears, and a hearing loss in their better ear.

### How it works

This option sends sound from a microphone, placed in your ear with little to no hearing, to a hearing aid receiver placed in your better-hearing ear. This allows you to hear from both sides. Your hearing care professional can help you determine what solution will best fit your needs and lifestyle.

# Say hello to the first rechargeable hearing aid for single- sided loss



Our new best-in-class A4 iQ CROS System rechargeable hearing aids last 20 percent longer than other systems, are 30 percent smaller and are 100 percent easy to use. Plus, with the industry's longest lasting charge, your hearing aids will always be ready to go.




*Our best hearing technology just got better.*

- **Immersive sound** for true listening enjoyment
- **Longest lasting charge** delivers over 30 hours of superior hearing
- **All-in-one charger** holds enough charge to provide portable charging without plugging into the wall
- **Smallest rechargeable** hearing aid available today

# Say hello to convenient connectivity

With more and more wireless products on the market, making connections between devices can be challenging. To navigate listening environments with high wireless traffic, our 900sync™ technology provides a powerful solution. Its high bit-rate streaming allows for ear-to-ear phone streaming and wireless capabilities for single-sided hearing loss.



## SURFLINK™

*Our Audibel wireless hearing aids wouldn't be complete without our popular SurfLink accessories.*



Our compact mobile phone transmitter, **SurfLink™ Mobile 2**, is an assistive listening device, media streamer and hearing aid remote all rolled into one.



**SurfLink Mini Mobile** is a small and simple mobile phone streaming solution. With this device, you can connect any Bluetooth-ready device to your A4 iQ hearing aids.



**SurfLink Media 2** is a set-and-forget media streaming solution that connects to TVs, MP3 players and more, to wirelessly stream audio directly to your hearing aids.



Our **SurfLink Remote** can be used to adjust volume, change memory modes, go in and out of streaming mode and much more — all with the push of a button.



Our discreet, lightweight **SurfLink Remote Microphone 2** is worn by a conversation partner to aid patients in one-on-one conversations in difficult listening environments.

# Say hello to the very best

## You want

*To hear comfortably  
in challenging  
environments*

*Rechargeable  
hearing solutions*

*To easily adapt  
to new sounds*

*Better wireless  
streaming*

*Natural one-on-one  
conversations*

*A comfortable,  
personalised listening  
experience*

*No buzzing or  
whistling*

*Durable,  
dependable  
hearing aids*

*Music the way it was  
meant to be heard*

*Customisable tinnitus relief*

## A4 iQ delivers

**Acuity™ Immersion Directionality**, designed to restore natural benefits provided by the ear and improve speech audibility in difficult listening situations.

**NEW** New best-in-class rechargeable hearing aids ensure you are always ready to go.\*

This feature speeds your transition to hearing aids by gradually adjusting settings to allow your brain to get used to new sounds.

Technology that provides consistent wireless performance while streaming mobile phones, TV, music and other media with **SurfLink Mobile 2**, **SurfLink Mini Mobile** and other **SurfLink** accessories.

The **SurfLink Remote Microphone** is a discreet, lightweight microphone worn by a conversation partner to aid patients in one-on-one conversations in difficult listening environments.

Capabilities that speed your transition to hearing aids with gradual adjustments, distortion-free listening in noise, and ultimate clarity for soft sounds.

Audibel's best-in-class feedback cancellation system providing feedback-free, comfortable listening all-day long.

**Surface™ NanoShield**, our pioneering water, wax and moisture repellent system, to protect and ensure durability and dependability.

A design that allows you to experience music in a whole new way. Now you can hear every note the way the artist intended — with pure, refined sound quality.

Advanced **Multiflex Tinnitus Technology** integrated into A4 iQ to bring relief to those who suffer from ringing in the ears.

# Say hello to better hearing everywhere you are

*From the beginning,  
our hearing solutions  
have enabled people to  
hear better in challenging  
environments and  
live a better life.*

You want to hear all the subtleties of life: the inflection in the voice of a spouse or the nuances in the notes of a favourite song. Whether it's one-on-one or in the most crowded places, you'll hear the world around you.

What you need will depend on how busy and active you are, and what types of environments you frequent.



*What is  
your hearing  
lifestyle?*

Our lineup of products enhances many different listening environments. Talk to your hearing professional about which options fit your lifestyle best.

*Compare the listening  
environments on the  
chart inside.*

# Say hello to better hearing everywhere you are

*From the beginning, our hearing solutions have enabled people to hear better in challenging environments and live a better life.*

You want to hear all the subtleties of life: the inflection in the voice of a spouse or the nuances in the notes of a favourite song. Whether it's one-on-one or in the most crowded places, you'll hear the world around you.

What you need will depend on how busy and active you are, and what types of environments you frequent.



## Listening environments

Listening Environments

Technology Levels



# Say hello to your unique style

## *Styles and colours*



AUDIBEL <sup>A</sup><sub>4</sub> iQ

Pick your favourite colour among our new RIC and BTE styles, both with tinnitus and wireless CROS and BiCROS solutions.

Consult your hearing care professional to see if rechargeable hearing aids are right for you.

### Colour guide



Black



Slate



Sterling



Bright White with Sterling



Espresso



Bronze



Champagne



*A4 iQ is available in a variety of styles and colours for your ideal fit.*

# Say hello to surround sound.

---

*Audibel delivers the best possible sound to your better ear so you don't have to worry about repositioning yourself for conversations or missing important sounds.*

Single-sided hearing loss doesn't have to get in the way of your lifestyle or your activities. With the A4 iQ CROS System, you can experience sound from all sides — never missing a single moment that brings you joy. At Audibel, we believe that to be fit to hear is to be fit for life.



AUDIBEL **A4** iQ



# Made so the world may hear



*Our Founder and CEO,  
William F. Austin, has  
always held a strong belief  
that the gift of better  
hearing results in a fuller  
and richer life.*

We prove our dedication to this philosophy daily, through our support of Starkey Hearing Foundation.

The gift of hearing opens up new worlds of opportunity. It connects individuals to life and helps them realise they can accomplish more than they ever thought possible. Starkey Hearing Foundation uses hearing as a vehicle to change lives around the world.



*Made to help millions.*



**100+**  
COUNTRIES



**MORE THAN  
1 MILLION  
PEOPLE IN NEED  
SERVED**



**COUNTLESS  
LIVES CHANGED**

Spread hope from ear to ear. To learn more or give the gift of hearing, visit [StarkeyHearingFoundation.org](http://StarkeyHearingFoundation.org).



[www.audibel.com](http://www.audibel.com)

---